

Pe Rabbit Anti-Human Trop-2 Monoclonal Antibody For Flow Cytometry - P09758

Item Number: CM0008612



Introduction

PE-conjugated rabbit monoclonal antibody targeting human TROP-2 for flow cytometry applications. Recognizes the tumor-associated antigen TROP-2, supporting immuno-oncology research and diagnostic development. - P09758

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human TROP-2 mAb
Remarks/Alias	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1
Species	Human
Gene ID	4070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-274 of human TROP-2 (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09758
Protein Weight	36kDa
Shipping	Ice bag