

Pe Rabbit Anti-Human/Monkey Cd10 Mab For Fc - P08473

Item Number: CM0010951



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD10 (Neprilysin) for flow cytometry. Cross-reacts with human and cynomolgus monkey. Validated for immunological and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD10 mAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43
Species	Human
GeneID (Human)	4311
GeneID	4311
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of human CD10 (NP_000893.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08473
Protein Weight	86kDa
Shipping	Ice bag