

Pe Rabbit Anti-Human/Monkey Cd127/IL-7R Monoclonal Antibody For Flow Cytometry - P16871

Item Number: CM0009431



Introduction

PE-conjugated rabbit monoclonal antibody against human/monkey CD127/IL-7R for flow cytometry. Recognizes the interleukin-7 receptor alpha chain, cross-reacts with human and cynomolgus monkey. Suitable for immunology research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD127/IL-7R mAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha
Species	Human
GeneID (Human)	3575
GeneID	3575
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-236 human CD127/IL-7Rα (NP_002176.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag