


Pe Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Flow Cytometry - P08571

Item Number: CM0008826



Introduction

PE-conjugated rabbit monoclonal antibody against human and cynomolgus monkey CD14, validated for flow cytometry. CD14 is a coreceptor for bacterial LPS, mediating innate immunity. Ideal for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID (Human)	929
GeneID	929
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-352 of human CD14 (NP_000582.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag