

Pe Rabbit Anti-Human/Monkey Cd16 Monoclonal Antibody For Fc - P08637/O75015

Item Number: CM0007806



Introduction

PE-conjugated rabbit monoclonal antibody targeting human/monkey CD16 (FCGR3A/FCGR3B). Validated for flow cytometry (FC), cross-reacts with human and cynomolgus. Ideal for immune cell profiling and ADCC studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD16 mAb
Remarks/Alias	CD16; FCG3; CD16A; FCGR3; IGFR3; IMD20; FCR-10; FCR11; CD16-II; FCGR11; FCR11A; FcGR11A; CD16b
Species	Human
GeneID (Human)	2214/2215
GeneID	2214/2215
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08637/O75015
Protein Weight	29KDa/26KDa
Shipping	Ice bag