

Pe Rabbit Anti-Human/Monkey Cd18 Monoclonal Antibody For Flow Cytometry - P05107

Item Number: CM0008649



Introduction

PE-conjugated rabbit monoclonal antibody against human CD18 (Integrin beta-2) for flow cytometry, cross-reactive with cynomolgus monkey. Ideal for studying leukocyte adhesion receptors and immune cell trafficking.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD18 mAb
Remarks/Alias	LAD; CD18; MF17; MF17; LCAMB; LFA-1; MAC-1
Species	Human
GeneID	3689
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-769 of human CD18 (NP_000202.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P05107
Protein Weight	85kDa
Shipping	Ice bag