

Pe Rabbit Anti-Human/Monkey Cd19 Monoclonal Antibody - Cd19

Item Number: CM0005012

Introduction

PE-conjugated rabbit monoclonal antibody targeting human and monkey CD19, validated for flow cytometry. Ideal for immunophenotyping and B-cell research in diagnostic and life science applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD19 mAb
Remarks/Alias	CD19
Species	Monkey
Gene ID (Human)	930
Gene ID (Monkey)	102145514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19 (XP_005591597.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Shipping	Ice bag