

Pe Rabbit Anti-Human/Monkey Cd20 Monoclonal Antibody For Fc - P11836

Item Number: CM0011218



Introduction

PE-conjugated rabbit monoclonal antibody against human/monkey CD20 (MS4A1), validated for flow cytometry. Specific for B-lymphocytes, with cross-reactivity to cynomolgus. Research use for B-cell studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD20 mAb
Remarks/Alias	MS4A1; B1; Bp35; CD20; CVID5; LEU-16; MS4A2; S7; B-lymphocyte antigen CD20
Species	Human
GeneID (Human)	931
GeneID	931
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 (NP_068769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P11836
Protein Weight	14kDa/33kDa
Shipping	Ice bag