

Pe Rabbit Anti-Human/Monkey Cd25 Monoclonal Antibody For Flow Cytometry - P01589

Item Number: CM0008858



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD25 (IL2RA, p55). Validated for flow cytometry on human and cynomolgus samples. Ideal for Treg and immune tolerance research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD25 mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
Gene ID (Human)	3559
Gene ID	3559
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag