

Pe Rabbit Anti-Human/Monkey Cd3 Mab For Fc - P07766

Item Number: CM0009956



Introduction

PE-conjugated rabbit monoclonal antibody targeting human/monkey CD3 for flow cytometry. Cross-reacts with cynomolgus CD3. Part of the TCR complex, essential for T-cell activation and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
Gene ID (Human)	916
Gene ID	916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag