

Pe Rabbit Anti-Human/Monkey Cd4 Monoclonal Antibody For Flow Cytometry - P01730

Item Number: CM0010170



Introduction

PE-conjugated rabbit monoclonal antibody targeting human/monkey CD4 for flow cytometry (FC). Validated for human and cynomolgus reactivity. 51kDa. Explore T-cell coreceptor function and HIV receptor studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4, IMD79, Leu-3, OKT4D, CD4mut
Species	Human
GeneID	920
Immunogen	Recombinant protein of human CD4
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag