

Pe Rabbit Anti-Human/Monkey Cd64 Monoclonal Antibody For Flow Cytometry - P12314

Item Number: CM0008729



Introduction

High affinity PE-conjugated rabbit monoclonal antibody against human and monkey CD64 (FcγRI), validated for flow cytometry. Ideal for studying monocyte/macrophage-specific IgG receptor functions in innate and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD64 mAb
Remarks/Alias	CD64; FCG1; FCRI; CD64A; FCGR1; IGFR1; FcgammaRI
Species	Human
Gene ID	2209
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P12314
Protein Weight	43kDa
Shipping	Ice bag