

Pe Rabbit Anti-Human/Monkey Cd69 Monoclonal Antibody For Fc - Q07108

Item Number: CM0010301



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD69, a key early activation antigen on immune cells. Validated for flow cytometry (FC), cross-reacts with cynomolgus monkey. Ideal for studying immune cell activation and regulatory T-cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey CD69 mAb
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26
Species	Human
GeneID (Human)	969
GeneID	969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-199 of human CD69 (NP_001772.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag