

Pe Rabbit Anti-Human/Monkey Hla-Dr Monoclonal Antibody For Flow Cytometry - P01903

Item Number: CM0010393



Introduction

PE-conjugated rabbit monoclonal antibody against human/monkey HLA-DR, validated for flow cytometry. Recognizes the alpha chain of MHC class II, expressed on antigen-presenting cells. Suitable for immunology and non-human primate research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human
GeneID (human)	3122
GeneID	3122
Immunogen	Recombinant protein of human HLA-DR.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag