


Pe Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Fc - P01880

Item Number: CM0012111



Introduction

PE-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey IgD (constant region, amino acids 99-384). Validated for flow cytometry, ideal for B-cell immunophenotyping and immunoglobulin receptor studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey IgD mAb
Remarks/Alias	IGHD
Species	Human
GeneID (Human)	3495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01880
Protein Weight	47kDa
Shipping	Ice bag