

Pe Rabbit Anti-Human/Monkey IgG Mab - Aaq57561.1

Item Number: CM0009862

Introduction

PE-conjugated rabbit monoclonal antibody against human and monkey IgG, suitable for intracellular flow cytometry. Immunogen: recombinant protein covering amino acids 98-322 of monkey IgG (AAQ57561.1). Cross-reacts with human and cynomolgus IgG.

[Learn More](#)

Parameter	Details
Product Name	PE Rabbit anti-Human/Monkey IgG mAb
Species	Monkey
Immunogen	Recombinant fusion protein containing amino acids 98-322 of monkey IgG (AAQ57561.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Cynomolgus
Shipping	Ice bag