


Pe Rabbit Anti-Human/Monkey Igm Mab For Fc - P01871

Item Number: CM0008776



Introduction

PE-conjugated rabbit monoclonal antibody targeting human/monkey IgM constant region, suitable for flow cytometry. Research tool for B cell immunology and humoral immunity studies. Recognizes IGHM (P01871).

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Monkey IgM mAb
Remarks/Alias	MU; VH; AGM1
Species	Human
GeneID	3507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 218-453 of human IgM (P01871).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01871
Protein Weight	49kDa
Shipping	Ice bag