

## Pe Rabbit Anti-Human/Monkey Ki67 Monoclonal Antibody For Intracellular Flow Cytometry - P46013

Item Number: CM0011200



### Introduction

PE-conjugated rabbit monoclonal antibody targeting human/monkey Ki67 (P46013), validated for intracellular flow cytometry. Essential tool for cell proliferation assays, cell cycle analysis, and oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                     |
|------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human/Monkey Ki67 mAb                                                                      |
| Remarks/Alias    | KiA, MIB-, MIB-1, PPP1R105                                                                                |
| Species          | Human                                                                                                     |
| GeneID (human)   | 4288                                                                                                      |
| GeneID           | 4288                                                                                                      |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 2271-2370 of human Ki67 (NP_002408.3). |
| Source           | Rabbit                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                     |
| Application      | FC (intra)                                                                                                |
| Cross Reactivity | Human, Cynomolgus                                                                                         |
| SWISS            | P46013                                                                                                    |
| Protein Weight   | 319kDa/359kDa                                                                                             |
| Shipping         | Ice bag                                                                                                   |