

Pe Rabbit Anti-Human/Mouse Integrin Alpha 5 (Itga5/Cd49E) Monoclonal Antibody For Flow Cytometry - P08648

Item Number: CM0012693



Introduction

PE-conjugated rabbit monoclonal antibody targeting human/mouse integrin alpha 5 (ITGA5/CD49e, VLA-5). Validated for flow cytometry. Ideal for cell adhesion, fibronectin receptor, and virus entry research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Mouse Integrin alpha 5 (ITGA5/CD49e) mAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A
Species	Human
Gene ID	3678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-995 of human Integrin alpha 5 (ITGA5/CD49e) (NP_002196.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag