


Pe Rabbit Anti-Human/Mouse Nanog Mab - Q80Z64

Item Number: CM0009320



Introduction

PE-conjugated rabbit monoclonal antibody targeting human and mouse Nanog, a homeobox transcription factor critical for pluripotency. Designed for intracellular flow cytometry (FC) in stem cell research and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Mouse Nanog mAb
Remarks/Alias	ENK; Stm1; ecat4; 2410002E02Rik
Species	Mouse
Gene ID (Human)	79923
Gene ID	71950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-305 of mouse Nanog (NP_082292.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Mouse
SWISS	Q80Z64
Protein Weight	34kDa
Shipping	Ice bag