

Pe Rabbit Anti-Human/Mouse Pdgfr α /Cd140A Monoclonal Antibody For Fc - P16234

Item Number: CM0010141



Introduction

PE-conjugated rabbit monoclonal antibody targeting human and mouse PDGFR α /CD140a (SWISS: P16234), validated for flow cytometry. Ideal for cell proliferation, embryonic development, and oncology research.

[Learn More](#)

Parameter	Description
Product Name	PE Rabbit anti-Human/Mouse PDGFR α /CD140a mAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2
Species	Human
Gene ID (Human)	5156
Gene ID	5156
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-524 of human PDGFR α /CD140a (NP_006197.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag