


Pe Rabbit Anti-Human/Mouse Phospho-Zap70-Y319+Syk-Y352 Monoclonal Antibody For Intracellular Flow Cytometry - P43403

Item Number: CM0000922



Introduction

PE-conjugated rabbit monoclonal antibody targeting phosphorylated ZAP70 (Y319) and Syk (Y352), validated for intracellular flow cytometry in human and mouse. Ideal for T-cell and B-cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Mouse Phospho-Zap70-Y319+Syk-Y352 mAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70
Species	Human, Mouse
Gene ID (Human)	7535
Gene ID	7535
Immunogen	A synthetic phosphorylated peptide around Y319 of human ZAP70 (NP_001070.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	FC (intra)
Cross Reactivity	Human, Mouse
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag