

Pe Rabbit Anti-Human/Mouse/Monkey Foxp3 Mab - Q9Bzs1

Item Number: CM0010395



Introduction

PE-conjugated rabbit monoclonal antibody targeting FOXP3 for intracellular flow cytometry. Validated in human, mouse, cynomolgus. Detects a 47kDa forkhead transcription factor essential for Treg-mediated immune regulation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human/Mouse/Monkey FOXP3 mAb
Remarks/Alias	JM2; AIID; IPEX; PIDX; XPID; DIETER
Species	Human
GeneID	50943
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXP3 (NP_054728.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human, Mouse, Cynomolgus
SWISS	Q9BZS1
Protein Weight	47kDa
Shipping	Ice bag