

Pe Rabbit Anti-Mouse Ccr6/Cd196 Monoclonal Antibody For Fc - O54689

Item Number: CM0005982

Introduction

PE Rabbit anti-Mouse CCR6/CD196 monoclonal antibody validated for flow cytometry (FC). Specific for mouse CCR6, a chemokine receptor critical for immune cell migration and mucosal immunity. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CCR6/CD196 mAb
Remarks/Alias	CCR-6; KY411; Cmkbr6; CC-CKR-6
Species	Mouse
GenelD (Human)	1235
GenelD	12458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-38 of mouse CCR6/CD196 (NP_033965.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O54689
Protein Weight	42kDa
Shipping	Ice bag