

Pe Rabbit Anti-Mouse Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P05532

Item Number: CM0010059



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD117/c-Kit. Validated for flow cytometry. Recognizes the mast/stem cell growth factor receptor Kit (P05532). Ideal for hematopoietic stem cell and mast cell research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD117/c-Kit mAb
Remarks/Alias	W; Bs; Fdc; Ssm; SC01; SC05; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3
Species	mouse
GeneID (human)	3815
GeneID	16590
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117/c-Kit (NP_001116205.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	22kDa/109kDa
Shipping	Ice bag