

Pe Rabbit Anti-Mouse Cd11A Monoclonal Antibody For Fc - P24063

Item Number: CM0005842

Introduction

PE Rabbit anti-Mouse CD11a monoclonal antibody validated for flow cytometry (FC) with mouse samples; targets integrin alpha-L involved in leukocyte adhesion and immune response - suitable for immunological research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD11a mAb
Remarks/Alias	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Species	Mouse
Gene ID (Human)	3683
Gene ID	16408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 149-340 of mouse CD11a (NP_001240801.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24063
Protein Weight	128kDa
Shipping	Ice bag