


Pe Rabbit Anti-Mouse Cd11B Monoclonal Antibody For Fc - P05555

Item Number: CM0008321



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD11b (ITGAM, MAC-1), validated for flow cytometry. Reacts with mouse monocytes, granulocytes, and microglia. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik
Species	Mouse
Gene ID (human)	3684
Gene ID	16409
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	115kDa/127kDa
Shipping	Ice bag