


Pe Rabbit Anti-Mouse Cd122/IL-2R β Monoclonal Antibody For Flow Cytometry - P16297

Item Number: CM0009978



Introduction

PE-labeled rabbit anti-mouse CD122/IL-2R β monoclonal antibody for flow cytometry. Detects the IL-2 receptor beta subunit, essential for IL-2 signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD122/IL-2R β mAb
Remarks/Alias	p70; CD122; IL15Rbeta; IL-2Rbeta; IL-15Rbeta; IL-2/15Rbeta
Species	Mouse
Gene ID (Human)	3560
Gene ID	16185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-235 of mouse CD122/IL-2R β (NP_032394.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16297
Protein Weight	60 kDa
Shipping	Ice bag