

Pe Rabbit Anti-Mouse Cd124/Il-4Ra Monoclonal Antibody For Fc - P16382

Item Number: CM0007765

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD124/IL-4R α , validated for flow cytometry (FC). Suitable for studying IL-4/IL-13 signaling in Th2 differentiation and allergic inflammation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD124/IL-4Ra mAb
Remarks/Alias	Il4r; CD124
Species	Mouse
Gene ID (Human)	3566
Gene ID (Mouse)	16190
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-233 of mouse CD124/IL-4Ra (NP_001008700.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16382
Protein Weight	87kDa
Shipping	Ice bag