

Pe Rabbit Anti-Mouse Cd127/Il-7Rα Monoclonal Antibody For Flow Cytometry - P16872

Item Number: CM0011446



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD127/IL-7 receptor alpha, validated for flow cytometry. Ideal for mouse immune cell identification and characterization in research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD127/IL-7Rα mAb
Remarks/Alias	CD127; IL-7Ralpha
Species	Mouse
GeneID (Human)	3575
GeneID	16197
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of mouse CD127/IL-7Rα(NP_612564.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16872
Protein Weight	51kDa
Shipping	Ice bag