


Pe Rabbit Anti-Mouse Cd130/Gp130 Monoclonal Antibody For Flow Cytometry - Q00560

Item Number: CM0004738



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD130/gp130 (IL-6 receptor subunit beta). Validated for flow cytometry (FC) with mouse samples. Core component of the IL-6 signaling pathway for immune and inflammatory research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD130/gp130 mAb
Remarks/Alias	CD130; gp130; D13Ert699e; 5133400A03Rik
Species	Mouse
GeneID (Human)	3572
GeneID	16195
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-617 of mouse CD130/gp130 (NP_034690.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00560
Protein Weight	102kDa
Shipping	Ice bag