

Pe Rabbit Anti-Mouse Cd138/Syndecan-1 Monoclonal Antibody For Flow Cytometry - P18828

Item Number: CM0010115

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD138/Syndecan-1 (SWISS: P18828). Validated for flow cytometry (FC) applications, with no cross-reactivity to other species. Recombinant immunogen covering extracellular domain. Ideal for immunophenotyping and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD138/Syndecan-1 mAb
Remarks/Alias	Sstn; Synd; CD138; Synd1; syn-1; CD138/Syndecan-1
Species	Mouse
Gene ID (Human)	6382
Gene ID (Mouse)	20969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-255 of mouse CD138/Syndecan-1 (NP_035649.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18828
Protein Weight	33kDa
Shipping	Ice bag