


Pe Rabbit Anti-Mouse Cd14 Monoclonal Antibody For Flow Cytometry - P10810

Item Number: CM0009589



Introduction

PE Rabbit anti-Mouse CD14 monoclonal antibody for flow cytometry, targeting monocyte differentiation antigen CD14. Suitable for innate immune response, LPS signaling, and inflammation research. Mouse cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD14 mAb
Remarks/Alias	CD14
Species	Mouse
GeneID (Human)	929
GeneID	12475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-345 of mouse CD14 (NP_033971.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P10810
Protein Weight	39kDa
Shipping	Ice bag