

Pe Rabbit Anti-Mouse Cd144/Ve-Cadherin Monoclonal Antibody For Flow Cytometry - P55284

Item Number: CM0010533



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse VE-cadherin (CD144), clone 7B4. Validated for flow cytometry applications in mouse samples. Facilitates research into endothelial cell adhesion, vascular permeability, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD144/VE-Cadherin mAb
Remarks/Alias	7B4; Vec; VECD; Cd144; VEcad; VE-Cad
Species	Mouse
GeneID (Human)	1003
GeneID	12562
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 46-592 of mouse CD144/VE-Cadherin (NP_033998.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P55284
Protein Weight	88kDa
Shipping	Ice bag