

Pe Rabbit Anti-Mouse Cd150/Slam Monoclonal Antibody For Fc - Q9Qum4

Item Number: CM0011753

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD150/SLAM, validated for flow cytometry (FC). Ideal for investigating SLAM-mediated immune cell activation, differentiation, and signaling pathways.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Mouse CD150/SLAM mAb
Remarks/Alias	Slam; CD150; IPO-3; CDw150; ESTM51; 4933415F16
Species	Mouse
GeneID (Human)	6504
GeneID	27218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-242 of mouse CD150/SLAM (NP_038758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QUM4
Protein Weight	36kDa/38kDa
Shipping	Ice bag