

Pe Rabbit Anti-Mouse Cd152/Ctla-4 Monoclonal Antibody For Flow Cytometry - P09793

Item Number: CM0012289



Introduction

PE-conjugated rabbit monoclonal antibody against mouse CTLA-4 (CD152). Designed for flow cytometry, enabling analysis of this T-cell inhibitory receptor in lymphoid tissues and immune checkpoint studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD152/CTLA-4 mAb
Remarks/Alias	Ctla4, Cd152, CTLA4
Species	Mouse
Gene ID (Human)	1493
Gene ID	12477
Immunogen	Recombinant protein corresponding to amino acids 36-161 of mouse CD152 (NP_033973.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P09793
Protein Weight	24kDa
Shipping	Ice bag