

Pe Rabbit Anti-Mouse Cd155/Pvr Monoclonal Antibody For Flow Cytometry - Q91Wp1

Item Number: CM0010746

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD155/Pvr (PVS, Taa1, necl-5). Validated for flow cytometry without secondary reagent. Recognizes native receptor. Molecular weight ~45kDa. Cross-reacts with mouse. Ideal for immunophenotyping and cancer studies.

[Learn More](#)

Specification	Value
Product Name	PE Rabbit anti-Mouse CD155/Pvr mAb
Remarks/Alias	PVS; mE4; HVED; Taa1; CD155; Tage4; necl-5; D7Ert458e; 3830421F03Rik
Species	Mouse
GeneID (Human)	5817
GeneID	52118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-348 of mouse CD155/Pvr (NP_612564.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q91WP1
Protein Weight	45kDa
Shipping	Ice bag