

Pe Rabbit Anti-Mouse Cd161C/Nk1.1 Monoclonal Antibody For Fc - P27814

Item Number: CM0010244

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD161c/NK1.1 (clone NKR-P1C). Validated for flow cytometry, detects mouse NK cells. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD161c/NK1.1 mAb
Remarks/Alias	Nk1; Ly59; NK-1; CD161; Ly-59; Ly55c; NK1.1; NKR-P1; Nk1.2; Klrb1b; NK-RP1; Nk-1.2; Nkrp1b; Nkrp1c; ly-55c; NKR-P140; NKR-P1.9
Species	Mouse
Gene ID (Human)	3820
Gene ID	17059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-223 of mouse CD161/NK1.1 (NP_001153376.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P27814
Protein Weight	25kDa
Shipping	Ice bag