

Pe Rabbit Anti-Mouse Cd163 Mab For Flow Cytometry - Q2Vlh6

Item Number: CM0011744

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD163 (M130), suitable for flow cytometry (FC). Targets the scavenger receptor cysteine-rich type 1 protein M130 involved in hemoglobin clearance and anti-inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD163 mAb
Remarks/Alias	CD163v2; CD163v3; M130
Species	Mouse
Gene ID (Human)	9332
Gene ID (Mouse)	93671
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-260 of mouse CD163 (NP_444324.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q2VLH6
Protein Weight	121kDa/125kDa
Shipping	Ice bag