


Pe Rabbit Anti-Mouse Cd18 Monoclonal Antibody For Fc - P11835

Item Number: CM0003986



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD18 (integrin beta-2), validated for flow cytometry. Reacts specifically with mouse samples, 85 kDa. Ideal for leukocyte adhesion, immune cell migration, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD18 mAb
Remarks/Alias	2E6; LAD; Cd18; Lfa1; MF17; LCAMB
Species	Mouse
GeneID (Human)	3689
GeneID	16414
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-771 of mouse CD18 (NP_032430.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P11835
Protein Weight	85 kDa
Shipping	Ice bag