

Pe Rabbit Anti-Mouse Cd19 Mab For Fc - P25918

Item Number: CM0008034

Introduction

PE-conjugated Rabbit anti-Mouse CD19 monoclonal antibody validated for flow cytometry (FC). B-cell marker coreceptor; ideal for immunoassay development and B-cell research. Reacts with mouse. - P25918

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD19 mAb
Species	Mouse
Gene ID (Human)	930
Gene ID	12478
Immunogen	Recombinant protein of mouse CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag