


Pe Rabbit Anti-Mouse Cd1d Monoclonal Antibody For Flow Cytometry - P11609 / P11610

Item Number: CM0010023



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD1d (P11609/P11610), validated for flow cytometry to study lipid antigen presentation and NKT cell biology in immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD1d mAb
Remarks/Alias	Cd1a; Cd1d; CD1.1; Ly-38
Species	Mouse
GeneID (Human)	912
GeneID (Mouse)	12479/12480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-305 of mouse CD1d (NP_031665.2/NP_031666.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P11609/P11610
Protein Weight	39kDa/38kDa
Shipping	Ice bag