

Pe Rabbit Anti-Mouse Cd200/Ox2 Monoclonal Antibody For Fc - O54901

Item Number: CM0004288

Introduction

PE-labeled rabbit anti-mouse CD200/OX2 monoclonal antibody validated for flow cytometry. Detects OX-2 membrane glycoprotein involved in T-cell costimulation and myeloid cell regulation. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD200/OX2 mAb
Remarks/Alias	OX2; Mox2
Species	Mouse
Gene ID (Human)	4345
Gene ID	17470
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-232 of mouse CD200/OX2 (NP_034948.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O54901
Protein Weight	31kDa
Shipping	Ice bag