


## Pe Rabbit Anti-Mouse Cd204 Monoclonal Antibody For Fc - P30204

Item Number: CM0005748



### Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD204 (Macrophage scavenger receptor types I/II). Validated for flow cytometry (FC), suitable for studying innate immunity, atherosclerosis, and macrophage function. UniProt P30204.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Mouse CD204 mAb                                                                                        |
| Remarks/Alias    | MSR; SR-A; Scvr; MRS-A; MSR-A; SR-AI; SR-AII; Scara1                                                                  |
| Species          | Mouse                                                                                                                 |
| GeneID (human)   | 4481                                                                                                                  |
| GeneID           | 20288                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 350-458 of mouse CD204 (NP_001106797.1) |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | FC                                                                                                                    |
| Cross Reactivity | Mouse                                                                                                                 |
| SWISS            | P30204                                                                                                                |
| Protein Weight   | 50kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |