

Pe Rabbit Anti-Mouse Cd23 Monoclonal Antibody For Fc - P20693

Item Number: CM0012312

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD23 (Fcer2), designed for flow cytometry. Reliable reagent for studying IgE receptor-mediated B cell and macrophage functions. Monoclonal, species-specific.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD23 mAb
Remarks/Alias	CD23; Fce2; Fcer2; Ly-42
Species	Mouse
GeneID (human)	2208
GeneID	14128
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-331 of mouse CD23 (NP_038545.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P20693
Protein Weight	38kDa
Shipping	Ice bag