

Pe Rabbit Anti-Mouse Cd235A/Glycophorin A Mab For Fc - P14220

Item Number: CM0009015

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD235a (Glycophorin A), validated for flow cytometry. A major erythrocyte membrane protein critical for membrane stability and MN blood group antigens. SWISS: P14220.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD235a/Glycophorin A mAb
Remarks/Alias	GPA; CD235a
Species	Mouse
GeneID (Human)	2993
GeneID	14934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of mouse CD235a (NP_034499.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P14220
Protein Weight	18 kDa
Shipping	Ice bag