

Pe Rabbit Anti-Mouse Cd24 Monoclonal Antibody For Flow Cytometry - P24807

Item Number: CM0008288

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD24 (HSA, Ly-52, nectadrin), validated for flow cytometry. Suitable for immunophenotyping and cell differentiation studies in mouse lymphoid, myeloid, and erythroid cells.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
GeneID (human)	100133941
GeneID	12484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag