

Pe Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Fc - P01590

Item Number: CM0008492

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD25 (IL-2R α , p55), suitable for flow cytometry applications in studying regulatory T cells and immune tolerance mechanisms. This antibody is specific for the mouse IL-2 receptor alpha chain and is ideal for Treg identification and functional assays.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GeneID (Human)	3559
GeneID	16184
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	30kDa
Shipping	Ice bag