

Pe Rabbit Anti-Mouse Cd262/Trailr2 Monoclonal Antibody For Fc - Q9Qzm4

Item Number: CM0004375

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD262/TRAILR2 (Death receptor 5) for flow cytometry. Suitable for apoptosis, death receptor signaling, and cancer research in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD262/TRAILR2 mAb
Remarks/Alias	MK; DR5; Ly98; KILLER; TRICKB; TRAILR2; TRICK2A; TRICK2B
Species	Mouse
Gene ID (Human)	8795
Gene ID	21933
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 53-170 of mouse CD262/TRAILR2 (NP_064671.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QZM4
Protein Weight	42kDa
Shipping	Ice bag