

Pe Rabbit Anti-Mouse Cd267/Taci Monoclonal Antibody For Flow Cytometry - Q9Et35

Item Number: CM0009800

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD267/TACI (Tnfrsf13b) for flow cytometry applications. Detects a 27 kDa receptor mediating B- and T-cell stimulation. Specific for mouse.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD267/TACI mAb
Remarks/Alias	Taci, 1200009E08Rik
Species	Mouse
GeneID (Human)	23495
GeneID (Mouse)	57916
Immunogen	Recombinant fusion protein containing amino acids 5-129 of mouse CD267/TACI (NP_067324.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Mouse
SWISS	Q9ET35
Protein Weight	27kDa
Shipping	Ice bag